

Schedule Copy					
Item No.	Description	Quantity & Packing	Sample Required		
1	2	3	4		
1	Polyster base round	150	3		
	corner X-Ray Films	Pkt.	Pkt.		
	(6 1/2" X 8 1/2")	1x 50	1 x 50		
	Green sensitive, blue	Films	Films		
	base with Base Fog				
	less than or equal to				
	0.23 & Dmax not less				
	than 3.10 Interleaved.				
	The density values should				
	be met at developer				
	temperature 32°C				
	to 35°C				
2	Polyster base round	9000	3		
	corner X-Ray Films	Pkt.	Pkt.		
	(8" X 10")	1x 50	1 x 50		
	Green sensitive, blue	Films	Films		

	base with Base Fog				
	less than or equal to				
	0.23 & Dmax not less				
	than 3.10 Interleaved				
	The density values should				
	be met at developer				
	temperature 32 ⁰ C				
	to 35 ⁰ C				
3	Polyster base round	17524	3		
	corner X-Ray Films	Pkt.	Pkt.		
	(10" X 12")	1x 50	1 x 50		
	Green sensitive, blue	Films	Films		
	base with Base Fog				
	less than or equal to				
	0.23 & Dmax not less				
	than 3.10 Interleaved				
	The density values should				
	be met at developer				
	temperature 32 ⁰ C				

	to 35°C				
4	Polyster base round	27855	3		
	corner X-Ray Films	Pkt.	Pkt.		
	(12" X 15")	1x 50	1 x 50		
	Green sensitive, blue	Films	Films		
	base with Base Fog				
	less than or equal to				
	0.23 & Dmax not less				
	than 3.10 Interleaved				
	The density values should				
	be met at developer				
	temperature 32°C				
	to 35°C				
	<u>X-RAY CHEMICALS</u>				
	<u>X-RAY EQUIPMENT</u>				
22	Wax-coated electro	7520	2		
	cardiographic paper	Rolls	Nos.		
	50 mm X 20 mts				
24	Lead Nos. 0 to 9	889	2		

	(1/2")	Sets	Sets		
25	Lead Letters	571	2		
	L & R 1/2"	Sets	Sets		
26	Lead Letters	453	2		
	A to Z 1/2"	Sets	Sets		
27	X-Ray Intensifying	217	2		
	Screens International	Pairs	Pairs		
	speed class 400 or				
	more Green emitting				
	(12" X 15")				
28	X-Ray Intensifying	205	2		
	Screens International	Pairs	Pairs		
	speed class 400 or				
	more Green emitting				
	(10" X 12")				
29	X-Ray Intensifying	206	2		
	Screens International	Pairs	Pairs		
	speed class 400 or				
	more Green emitting				
	(8" X 10")				

31	Contact Jelly for .	30445	2		
	Ultrasonography	Nos.	Nos.		
	which should be	(1 x 250			
	1. Non allergic, Non	gms.)			
	greezy, Non Toxic.				
	2. 100% Airless				
	Skin contact				
32	Papers for CST	1200	2		
	Computerized stress	Pkts.	Pkts.		
	test machine :	(300			
	1) B4 Universal Thermal/	Sheets/			
	Paper.	Pkt.)			
	2) Full page Red grid				
	3) 216 x 280 mm				
	(8.5" x 11")				
	4) Z-fold				
	5) Queue hole				
	6) 300 Sheets /Pkt.				
33	Compact disk (CD-R)	38830	2		
	atleast 700MB each	Nos.	Nos.		

	with separate cover				
	at least 48X write				
	speed-				
34	DVD-R + atleast	20970	2		
	4GB atleast 4X write	Nos.	Nos.		
37	Sterile Braided	300	2		
	Angiographic catheters	Nos.	Nos.		
	4F to accept 0.038"				
	guide wire. 100 cm.				
	long. "H 1"				
38	Sterile Braided	900	2		
	Angiographic catheters	Nos.	Nos.		
	4F to accept 0.038"				
	guide wire. 60 cm.				
	long. "C 1"				
39	Sterile Braided	200	2		
	Angiographic catheters	Nos.	Nos.		
	4F to accept 0.038"				
	guide wire.100 cm.				
	long."Sim 1"				

40	Sterile Braided	320	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guide wire.100 cm.				
	long."Sim 1"				
41	Sterile Braided	300	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guide wire. 60 cm.				
	long. "C 1"				
42	Sterile Braided	350	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guide wire. 60 cm.				
	long. "RDC"				
43	Sterile Braided	300	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guide wire. 100 cm.				
	long. "H 1"				

44	Sterile Braided	300	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guidewire. 100cm long.				
	"Straight" with multiple				
	side holes				
45	Sterile Braided	600	2		
	Angiographic catheters	Nos.	Nos.		
	5F to accept 0.038"				
	guide wire. 100cm long.				
	"Pigtail with multiple				
	side holes.				
46	Vascular Introducer	450	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 4F (French)				
47	Vascular Introducer	2250	2		

	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 5F (French)				
48	Vascular Introducer	5900	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 6F (French)				
49	Vascular Introducer	2150	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 7F (French)				

50	Vascular Introducer	520	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 8F (French)				
51	Vascular Introducer	320	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				
	Size 9F (French)				
52	Vascular Introducer	290	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.038"/0.035" J. tip				
	guide wire for adults				
	(10 to 12 cm. long)				

	Size 10F (French)				
	<u>INJECTION</u>				
59	Inj. Meglumine	44350	10		
	Diatrizoate Sodium	Amp/	Amp/		
	Diatrizoate Iodine	Vial	Vial		
	content 365-375 mg/ml	(1 x 20	(1 x 20		
	in 20 ml. Ampoule/Vial	ml)	ml)		
60	Inj. Meglumine	13150	10		
	Diatrizoate Iodine	Amp/	Amp/		
	content 290-310 mg/ml	Vial	Vial		
	in 20 ml. Ampoule/Vial	(1 x 20	(1 x 20		
		ml)	ml)		
61	Inj. Meglumine	58840	5		
	Diatrizoate Iodine	Vial	Vial		
	content 290-310 mg/ml	(1 X 50	(1 X 50		
	in 50 ml. Vial	ml)	ml)		
62	Inj. Non ionic contrast	19960	5		
	for intravascular use	Vial	Vial		
	Iodine concentration	(1 x 100	(1 X 100		
	290-310 mg/ml	ml)	ml)		

	100 ml Vial.				
63	Inj. Non ionic contrast	44560	5		
	for intravascular use	Vial	Vial		
	Iodine concentration	(1 x 100	(1 X 100		
	350-370 mg/ml	ml)	ml)		
	100 ml Vial.				
64	Inj. Meglumine	31000	5		
	Diatrizoate Sodium	Vial	Vial		
	Diatrizoate Iodine	(1 x 50	(1 X 50		
	content 365-375 mg/ml	ml)	ml)		
	in 50 ml. Vial.				
65	Inj. Meglumine	15000	5		
	Diatrizoate Sodium	Vial	Vial		
	Diatrizoate Iodine	(1 x 50	(1 X 50		
	content 320-335 mg/ml	ml)	ml)		
	in 50 ml. Vial.				
75	Pigtail Angiography	850	2		
	catheter with multiple	Nos.	Nos.		
	side holes				
	Adult 6F				

76	Pigtail Angiography	275	2		
	catheter with multiple	Nos.	Nos.		
	side holes				
	Adult 7F				
77	Pigtail Angiography	450	2		
	catheter with multiple	Nos.	Nos.		
	side holes length 65 to				
	75 cm 5F				
78	Vascular Introducer	860	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.021" guide wire with				
	J. tip on one end &				
	straight soft tip on the				
	other end for Paediatrics				
	(7.5 cm long) Size 6F				
	(French)				
79	Vascular Introducer	1080	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				

	0.021" guide wire with				
	J. tip on one end &				
	straight soft tip on the				
	other end for Paediatrics				
	(7.5 cm long) Size 5F				
	(French)				
80	Vascular Introducer	620	2		
	sheath with valve,	Nos.	Nos.		
	sidearm, soft dilator &				
	0.021" guide wire with				
	J. tip on one end &				
	straight soft tip on the				
	other end for Paediatrics				
	(7.5 cm long) Size 4F				
	(French)				
81	Inj. Gadolinium based	6000	10		
	Contrast for MRI	Vial	Vial		
	for Paediatric use	1 x 10	1 x 10		
		ml.	ml.		
83	Inj. Gadolinium based	6000	10		

	Contrast for MRI	Vial	Vial		
	for general use	1 x 10	1 x 10		
		ml.	ml.		
84	Radiofrequency ablation	25	2		
	Catheter 4 mm tip small	Nos.	Nos.		
	curve (Green)				
85	Radiofrequency ablation	20	2		
	Catheter 4 mm tip	Nos.	Nos.		
	medium curve (Blue)				
86	Diagnostic nondeflatable	30	2		
	quadripolar catheter	Nos.	Nos.		
87	Diagnostic deflatable	15	2		
	decapolar catheter	Nos.	Nos.		
88	Diagnostic deflatable	10	2		
	duodecapolar catheter	Nos.	Nos.		
97	Chemical coated ECG Paper	4708	2		
	Single Channel	Rolls	Rolls		
	50 mm x 20 mtrs				
98	X-Ray Folders Laminated 250 gsm	190000	2		
	14 1/2" x 17 1/2"	Nos.	Nos.		

102	Dry View Kodak films cartridges	12000	2		
	of 125 films 14" x 17" for Kodak	Cartridges	Cartridges		
	Dry View Cameras				
103	Agfa Dry View Films Size 14" x 17"	500	3		
		Pkts.	Pkts.		
		(1 x 100			
		Films)			
104	Agfa Dry View Films Size 11" x 14"	500	3		
		Pkts.	Pkts.		
		(1 x 100			
		Films)			
105	Agfa Dry View Films Size 8" x 10"	800	3		
		Pkts.	Pkts.		
		(1 x 100			
		Films)			
106	ECG electrodes MR Compatible	9000	2		
		Nos.	Nos.		
122	Tissue Paper for cleaning	11000	2		
	Ultrasound Jelly.	Packs	Packs		
	(Size not less than 20 cm x 20 cm)	of 100			

		Papers			
125	Agfa Dry View Films	600	2		
	size 8" x 10"	Pkts.	Pkts.		
		1 x 100	1 X 100		
		Films	Films		
			Pkt.		
126	Konica Dry film for Dry	2000	2		
	Camera Size 8" x 10"	Pkts.	Pkts.		
		(Pkt. Of	Pkt. Of		
		125 films)	125 films		
127	Konica Dry film for Dry	2000	2		
	Camera Size 10" x 12"	Pkts.	Pkts.		
		(Pkt. Of	Pkt. Of		
		125 films)	125 films		
128	Konica Dry film for Dry	2000	2		
	Camera Size 14" x 17"	Pkts.	Pkts.		
		(Pkt. Of	Pkt. Of		
		125 films)	125 films		
Rate				Amount	Remarks
Basic	Taxes, Duties,	Total			Brand

	Octroi & Other	in Figures	in Words		Name Etc.
	Levies to be				
Rs. P.	shown separately	Rs. P.	Rs. P.		
6	7	8	9	10	11